

Aesculap[®] AdTec[®] comBi

Bipolar combination instrument for laparoscopy
4 in One... grasping – dissecting – coagulating – cutting.



Aesculap Endoscopic Technology

Aesculap® AdTec® comBi

4 in One... grasping – dissecting – coagulating – cutting.

4inOne

AdTec® comBi is a combination instrument for laparoscopic surgery. It combines a **bipolar grasping forceps** with **mechanical scissors**.


The integrated functions allow:

- precise dissection
- authentic grasping and manipulation of tissue
- a high level of dependability through bipolar coagulation
- effective mechanical cutting

Further advantages


- fine instrument tips for excellent overview
- colour coding for simple recognition of the active function (scissors or grasping forceps)
- "4 in 1" cost efficiency (grasping, dissecting, coagulating, cutting) – fewer instrument changes





Ergonomic handle:

A trigger integrated into the handle allows fast and easy switching from the grasping forceps to the scissors function.

			Cat. no.
<p>Bipolar combination instrument for single use, pack containing 6 pieces</p>	<p>5 mm/310 mm</p>	 <p>AESCULAP</p>	<p>PO800SU</p>
	<p>1:1</p>		
<p>The clinical outcomes demonstrate: special procedures in general surgery and gynaecology can be performed more accurately, faster and more cost effectively.</p>			

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany
Phone +49 7461 95-0 | Fax +49 7461 95-26 00 | www.aesculap.com

Aesculap – a B. Braun company

The main product trademark 'Aesculap' and the product trademark 'AdTec' are registered trademarks of Aesculap AG.

Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.

Brochure No. C24002

0612/2.0/4