

Bone Wax

Hemostat for Bone Bleeding



Bone Wax

Hemostat for Bone Bleeding


Sterile Mixture of Natural Beeswax and Vaseline

- Mechanical haemostasis in bones (1-5)
- Haemostasis achieved by tamponade effect (6-7)
- Achieve mechanical occlusion of intraosseous vessels (2)
- Soft, easy to shape and apply "only hand-warm" (2)



Indication

Bone Wax is indicated for the control of bleeding from bone surfaces during surgical procedures that involve bone cutting (8).



Ordering Information

Description	Code
24 Bars set	102 9754

- 1) Alhan C, Arıtırık C, Senay S, Okten M, Güllü AU, Kilic L et al. Use of bone wax is related to increased postoperative sternal dehiscence. *Kardiochir Torakochirurgia Pol.* 2014;11(4):385–90.
- 2) Zhou H, Ge J, Bai Y, Liang C, Yang L. Translation of bone wax and its substitutes: History, clinical status and future directions. *J Orthop Translat.* 2019;17:64–72.
- 3) Park JS, Kong D-S, Lee J-A, Park K. Intraoperative management to prevent cerebrospinal fluid leakage after microvascular decompression: dural closure with a "plugging muscle" method. *Neurosurg Rev.* 2007;30(2):139–42.
- 4) Serraino GF, Dominijanni A, Jiritano F, Rossi M, Cuda A, Caroleo S et al. Platelet-rich plasma inside the sternotomy wound reduces the incidence of sternal wound infections. *Int Wound J.* 2015;12(3):260–4.
- 5) Ma Q, Sun C-J, Wu S, Cai X. Comparison of Blood Loss between Open-Box Prosthesis and Closed-Box Prosthesis after Primary Total Knee Arthroplasty. *Orthop Surg.* 2021;13(3):768–77.
- 6) Schonauer C, Mastantuoni C, Somma T, Falco R de, Cappabianca P, Tessitore E. Topical hemostatic agents in neurosurgery, a comprehensive review: 15 years update. *Neurosurg Rev.* 2022;45(2):1217–32.
- 7) Achneck HE, Sileshi B, Jamiolkowski RM, Albala DM, Shapiro ML, Lawson JH. A comprehensive review of topical hemostatic agents: efficacy and recommendations for use. *Ann Surg.* 2010;251(2):217–28.
- (8) Instructions for use.



B. Braun Surgical, S.A. | Carretera de Terrassa, 121 | 08191 Rubí | Spain | Phone +34 93 586 6200 | www.bbraun.com